


DiaQ[®]

Blood Tubing Set



PRE PUMP

POST PUMP

Haemodialysis Blood tubing set has been specially designed for the patient to connect external system which extracts patient's blood to the dialyzer and then returns patient's blood from the dialyzer. Haemodialysis Blood Lines specially used for continuous blood supplies during haemodialysis process.

Consists of 2 parts: Arterial and Venous line which are used during dialysis with attached fistula and dialyzer.

FEATURES

- Highly kink resistant tubing
- Single-use, EO Sterilized & Individually Packed
- Large finger guard and latex free injection sites.
- Large and small clamps provide features for ideal ergonomics.
- To ensure secure lock of dialyser ports grip, connectors of wing shape provided.
- For secure and smooth connection to patient's vascular access, wing connectors of wing shape are provided.
- Drip chambers are available in different models and sizes, soft and rigid, straight and cone shaped as per the clients specifications.
- Medical grade PVC tube, ensure the biggest safety on human body and can easily equipped with different dialysis machines.
- Available in any configuration as per customer's requirement for various dialysis machines.



www.diacaresolutions.com



CORPORATE OFFICE :

Unit No. 408, 4th Floor, Eco Star, Vishweshwar Nagar, Off. Aarey Road, Goregaon (E), Mumbai - 400063.
Tel.: +91-22-45161475 | Email : enquiry@diacaresolutions.com

BRANCH OFFICE :

Office No. 212, Surya Kiran Building, Plot No 19, Connaught Place, K. G. Marg., New Delhi - 110001.
Mobile : +91-8800908729 | Email : enquiry@diacaresolutions.com