



DiflowTM

AV Fistula Needle



A.V. fistula needles set are indicated for use in conjunction with a connector of Hemodialysis blood tubing set. AV Fistula Needles used to connect blood lines to the blood vessel through needle when dialysis is carried out via an internal fistula. Flexible butterfly wings are for proper fixation. Color coded wings for easy identification of needle sizes. Soft, frosted and kink resistant PVC tubing during dialysis and Haemodialysis Process.

FEATURES

- Oval-shaped large back-eye prevent turbulence and minimize tissue damage.
- Color-coded bevel index for clear indication of needle position
- Premium luer-lock connector with wing for safe connection to the bloodline system
- Color-coded wings and packaging clearly indicating the gauge size
- Soft tubing with anti-kinking properties
- Excellent needle quality for smooth treatments
- Silicone coating to prevent blood coagulation.
- Sterilized in EO gas.



www.diacaresolutions.com



CORPORATE OFFICE :

Unit No. 408, 4th Floor, Eco Star, Vishweshwar Nagar, Off. Aarey Road, Goregaon (E), Mumbai - 400063.
Tel.: +91-22-45161475 | Email : enquiry@diacaresolutions.com

BRANCH OFFICE :

Office No. 212, Surya Kiran Building, Plot No 19, Connaught Place, K. G. Marg., New Delhi - 110001.
Mobile : +91-8800908729 | Email : enquiry@diacaresolutions.com